




The plot shows two spectral traces. The top trace is blue and shows a broad, noisy peak. The bottom trace is orange and shows a similar broad peak with a distinct dip in the center, characteristic of a saturation dip in spectroscopy.

Saturation Spectroscopy


Hyperfine Structure in ^{85}Rb



selection rule for
hyperfine
structure
transitions:
 $\Delta F=0, \pm 1$


Natural linewidth and Doppler broadening

- natural linewidth determined by the finite lifetime of the excited state, Lorentzian profile
- thermal velocity distribution causes Doppler shifts in the transition frequencies (thermal Doppler broadening)




Saturation Spectroscopy Set-up

- probe beam is overlapped with a counterpropagating pump beam with a higher intensity
- on resonance, pump and probe beam excite the same velocity class of atoms
- due to the higher intensity atoms are more likely to absorb photons from the pump beam → absorption minimum in the probe beam
- cross-over peaks for frequencies between transitions




Saturation Spectroscopy Set-up

- probe beam is overlapped with a counterpropagating pump beam with a higher intensity
- on resonance, pump and probe beam excite the same velocity class of atoms
- due to the higher intensity atoms are more likely to absorb photons from the pump beam → absorption minimum in the probe beam
- cross-over peaks for frequencies between transitions



Results for ^{85}Rb , D_2 , $F=3$


signals with and without the pump beam



Lamb dips


difference of the signals

Hyperfine structure transitions




cross-over peaks

Results for ^{85}Rb , D_2 , $F=3$




Hyperfine structure transitions



cross-over peaks

Application

- laser frequency stabilisation mechanism using a feedback proportional to the derivative of the obtained absorption signal
- utilizes the fixed frequency of a hyperfine structure transition
- feedback to the laser is zero on the resonance frequency, sign and magnitude for off-resonance frequencies depend on the derivation to the ideal frequency



Summary

- hyperfine structure splitting occurs due to interactions of the electrons with the atomic nucleus
- saturation spectroscopy relies on a pump-probe-setup to overcome thermal Doppler-broadening and can be used to probe hyperfine structure transitions at room temperature
- saturation spectroscopy setup can be used to generate a feedback signal to stabilise a laser on the fixed frequency of a hyperfine structure transition

Thanks to:

Dr. Pearson


Dr. Teigelhöfer

Prof. Buchinger


Prof. Crawford



Laser Stabilisation Set-up



Generation of the feedback signal



Absorption signal with frequency modulation (a) below resonance, (b) at resonance, and (c) above resonance.

In each plot the top curve is the absorption signal, the bottom continuous curve is the modulation, and the broken curve is the product.

Source:

A. Kumarakrishnan, M. Weel. "Laser-frequency stabilization using a lock-in amplifier". Can. J. Phys. 80: 1449–1458 (2002).